

402

Nitrile Powder-Free Non-Sterile Examination Glove



PRODUCT DESCRIPTION

Material	Nitrile
Colour	White
Shape	Ambidextrous
Cuff	Beaded
External surface	Finger textured
Internal surface	
Internal finish	Powder-Free

PHYSICAL PROPERTIES

Thickness (single)	Finger	min 0.060
Typical average values (mm)	Palm	min 0.050
	Cuff	min 0.040
	Strength (average values)	Before ageing
Elongation at break (%)	min 500	min 400
Force at break (N)	min 6.0	min 6.0
Minimum length (mm)	240	

Shelf life	3 years
Storage instructions	Keep out of direct sunlight. Store in a cool and dry place. Keep away from sources of ozone or ignition.
Packaging	20x100 pieces / 2000 pieces per case

PRODUCT REFERENCES

Size / product code	X-Small* 242GK265	Small 242GK270
	Medium 242GK275	Large 242GK280
	X-Large* 242GK285	

* for special orders

MANUFACTURING AND SAFETY STANDARDS

AQL (Pin-Holes)	Waterleak tested Conform European Norm EN 455-1: Inspection Level I AQL 1.5
Protein content	Contains no latex proteins
Sterilization	Not applicable
CE-classification Compliances	Class I CE mark EN 455 parts 1, 2, 3 - EN 374 parts 2,4,5
Food contact	Suitable for food contact EC regulation 1935/2004

KEY FEATURES AND ADVANTAGES

- Allergy prevention
- Strong safety and chemical resistance
- Comfort and ease-of-use

